Laboratory Shaker (Horizontal Shaking)

MK161

Rotary, elliptical and reciprocate motion

- Compact, space saving design
- Changeable rotary, elliptical and reciprocate motion for mixing, extracting and stirring of samples
- Stable and high torque shaking power and speed with the DC brushless motor
- Shaking frequency and timer are dial setting and digital display
- Shake pause function, timer function and constant operation by one switch
- Selectable mixing, extracting and stirring patterns when used with different shaking stage and racks (optional item)
- Can be placed inside IN602CSW incubator for shaking incubation

Specifications

opecifications				
MK161				
Rotary, Elliptical and reciprocate (Manual operation)				
Rotary:30mm Reciprocate: 30mm				
20~200rpm				
Dial setting, Digital display				
Dial setting, Digital display / Digital 0.1min. (6sec.) to 99.9hr.				
Main unit : W300 x D254mm, Stage : W290 x D250mm				
W350 x D300 x H150mm				
~15kg				
AC115V 0.5A / AC220V 0.3A				


Operational accessories

Mounting stage



Capacity	Number of erlenmeyer flask clamp
100ml	10pcs
200ml	9pcs
300ml	5pcs
500ml	4pcs
1,000ml	2pcs
Product code	232061

Diagonal centrifugal tube holder



*Mounting stage sold seperately

*Glassware not included

Erlenmeyer flask holder clamp



	Product code	Capacity	No. of clamps	
	232063	200ml	9 pcs.	
	232062	100ml	10 pcs.	
	232064	300ml	5 pcs.	
	232065	500ml	4 pcs.	
	232066	1,000ml	2 pcs.	
*Mounting stage sold seperately				

Non-skid sheet



 Product code
 Dimension (wx Dxir)

 232084
 290×250×30

 *Mounting stage sold seperately




Example of using mounting stage and erlenmeyer flask holder clamps (optional)



Incubator IN602CW slide shaker stage (optional) and MK161 installation example

*Glassware not included.

Diagonal rack holder



232082

232007	TUUIIII	5 pcs
232068	200ml	2 pcs
232069	300ml	2 pcs

Mounting stage sold seperately

Single spring shaking rack



 Dimension (W×D×H)

 290×250×66mm

 Number of test tube:

 ø16mm test tube×64 (45°inclination)

 Number of erlenmeyer flask

 50ml×20pcs, 100ml×10pcs, 200ml×9pcs,

 300ml×5pcs, 500ml×4pcs. 1000ml×2pcs

 Product code
 232050

*Mounting stage not necessary This can be set directly to the main unit.

Two layer spring shaking rack

ø18

Mounting stage sold seperately

20 pcs



 Dimension (W×D×H)

 290×250×110mm

 Number of test tube:

 ø16mm test tube×64 (45°inclination)

 Number of erlenmeyer flask

 50ml×20pcs, 100ml×10pcs, 200ml×9pcs,

 300ml×5pcs, 500ml×4pcs. 1000ml×2pcs

 Product code
 232056

*Mounting stage not necessary This can be set directly to the main unit.